



Megan P. Pinette<sup>1</sup>, Karen L. Troy<sup>1</sup>, Nguyen Nguyen<sup>2</sup>, Ricardo A. Battaglino<sup>2</sup>, Leslie R. Morse<sup>2</sup> <sup>1</sup>Department of Biomedical Engineering, Worcester Polytechnic Institute <sup>2</sup>Department of Physical Medicine and Rehabilitation, University of Minnesota

#### INTRODUCTION

- Spinal cord injuries (SCI) affect 28-55 million people in the United States [1].
- Following a SCI, rapid muscle atrophy leads to decreased metabolic rate, increased fat storage, and increased risk of obesity-related illnesses such as cardiovascular disease [2].
- In radiographs, adipose tissue has very low radio-density.
- American Spinal Injury Association (ASIA) classifies the level of injury based on an impairment scale: • Class A (Complete) through E (Normal)
- Robotic exoskeletons are becoming increasingly available in clinical settings. However, little is known about the musculoskeletal benefits of using an exoskeleton.
- Robotic exoskeletons rely on motors at the hip and knee, a computerized control system, and rechargeable batteries to allow for assisted walking and rehabilitation [4].
- Exoskeleton-assisted gait training can be used to deliver controlled repetitive training and reduce physical burden on physical therapists [3].

**Objective:** (1) To develop standardized methods to measure muscle radio-density and volume in CT scans. (2) To compare CT data between time points in an ongoing clinical trial.

### **METHODS**

#### **Subjects**

- 27 subjects with class A-C SCI (age: 33.7±8.8 years, height: 177.8 ± 10.2 cm, weight: 74.8 ± 14.1 kgs.) were enrolled in a 12 month clinical trial
- Subjects were randomly assigned to two groups:
  - Immediate gait training: subjects participated in exoskeleton assisted gait training for study months 0 to 6
  - Delayed gait training: subjects participated in exoskeleton assisted gait training for study months 6 to 12
- Gait training lasted for 1 hour, 3 days a week for 26 weeks

Immediate

0 months Baseline CT Scan - Begin Gait Training

6 months Mid-point CT Scan - End Gait Training; Usual Care

Delayed

Baseline CT Scan - Usual Care

Mid-point CT Scan - Begin Gait Training

# Exoskeleton-Assisted Walking Increases Muscle Volume in People with Spinal Cord Injury (SCI)

### METHODS CONT.



Figure 1. Indego Exoskeleton



#### **Quantitative Image Analysis**

- Mimics v. 18.0 (Materialise, Leuven, Belgium) was used to quantitatively analyze CT scans
- 15% regions of the thigh and 15% of the calf were cropped into individual sections on each leg, calculated based on subject height



**Figure 2**. 25-40% of the upper leg and 15-30% of the lower leg were cropped into individual sections. Analyzed regions in all three planes are shown in red.

- Muscles were isolated from skin, blood vessels, etc. using region growing, thresholding and morphology tools. A threshold of -49 to 169 Hounsfield units (HU) was used to delineate muscle tissue. Less radio-dense tissue (lower average HU) indicates fattier/lower quality muscle.

#### **Statistical Analysis**

- Pearson correlations were calculated for demographics and baseline values of muscle volume and radio-density Precision (%CV) calculated for all metrics
- Paired t-tests were used to compare pre/post measures of
- muscle volume, radio-density and volume\*radio-density

### RESULTS



Figure 3. Height (left) and mass (right) versus muscle volume.

- p=0.346)



Figure 4. Pre and post muscle radio-density (left) and volume (right). Left figure shows that muscle density did not significantly change during gait training, but right figure shows that muscle volume increased significantly by 15% (p=0.02). Red bars show the group average at each timepoint. Error bars represent standard deviation.

- quantity in people with SCI.

[1] J. Mcdonald and C. Sadowsky, "Spinal Cord Injury," The Lancet, vol. 359, no. 9304, pp. 417–425.[2] L. Giangregorio and N. McCartney, "Bone Loss and Muscle Atrophy in Spinal Cord Injury: Epidemiology, Fracture Prediction, and Rehabilitation Strategies," The Journal of Spinal Cord Medicine, vol. 29, no. 5, pp. 489–500.[3] S. Banala, S. Hun Kim, S. Agrawal, and J. Scholz, "Robot Assisted Gait Training With Active Leg Exoskeleton (ALEX)," IEEE *Transaction on Neural Systems and Rehabilitation Engineering*, vol. 17, no. 1.[4] C. Tefertiller, K. Hays, J. Jones, A. Jayaraman, C. Hartigan, T. Bushnik, and G. Forrest, "Initial Outcomes from a Multicenter Study Utilizing the Indego Powered Exoskeleton in Spinal Cord Injury.," Spinal Cord Injury Rehabilitation, vol. 24, no. 1, pp. 78–85, Nov. 2017.

# ACKNOWLEDGEMENTS

This work was funded by grant W81XWH-15-2-0078 (United States Department of Defense) to Leslie R. Morse. We thank Nguyen Nguyen for assisting with data collection, Leslie Morse for supervising the clinical trial, and all of the individual participants for volunteering their time.



### RESULTS

Precision for muscle volume and radio-density was 1.3-1.4%. Gait training was associated with significant increases in muscle radio-density\*volume (16% increase; p=0.035) These increases were primarily due to changes in volume (15% increase; p=0.02), and not radio-density (3% increase;

# DISCUSSION

This research demonstrates that passive motion generated during exoskeleton-assisted walking could improve muscle

Muscle mass is linked to improved cardiovascular fitness These results could be used to develop and evaluate improved rehabilitation interventions

# REFERENCES